



D E S K T O P

LAMINAR AIRFLOW

DESKTOP LAMINAR AIRFLOW LLFD-100 SERIES

Labocon Desktop Laminar Airflow LLFD-100 Series is specially designed for laboratory works and to provide sterility. This new series includes both horizontal and vertical laminar flow cabinet with a small volume. It is used in variety of applications in medical research laboratories, hospitals, manufacturing facilities and other research and production environments.

Application

Labocon Desktop Laminar Airflow is used in tissue culture, media preparation, medical and biological experiments and other research and production environments.

Specification

Model	LLFD-101	LLFD-102
Type	Horizontal Laminar Flow	Vertical Laminar Flow
Air Cleanliness	Class 100 @ $\geq 0.5\mu\text{m}$	
Average Air Flow Velocity	0.3-0.45m/s	
Fluorescent Lamp	8W	
Noise Level	$\leq 60\text{dB(A)}$	
Inner Dimension	600×495×500 mm	600×400×535 mm
Overall Dimension	650×535×860 mm	650×700×625 mm
Power	400W	
Power Supply	220 V, 60 Hz	
Catalog No.	9281124253	9281125253



LLFD-101