



Shaking Water Bath

Shaking Water Bath LSWB-100 Series



Labocon Shaking Water Bath LSWB-100 Series is equipped with microprocessor controller for temperature and speed. It offers accurate heating up processes and simultaneous gentle mixing or vigorously shaking to agitate samples.

Features:

- * Digital display.
- * Stainless steel internal chamber.
- * Easy to clean.
- * Equipped with clamp.



Application:

Shaking Water Bath LSWB-100 Series from Labocon is used in biology, biochemistry, industry and pharmacology applications. Their scope of applications also comprises in incubation, fermentation, homogenization, chemical reactions, enzyme studies and tissue studies.

Specification:

Model	LSWB-101	LSWB-102	LSWB-103
Capacity	50 Litres		
Method	Orbital	Reciprocal	Orbital and Reciprocal
Temperature Range	RT+5°C to 99°C		
Speed Range	0-220 rpm	0-300 rpm	0-220 rpm; 0-300 rpm
Orbit Diameter	30 mm	20 mm	30 mm; 20 mm
Timer Range	0-999 hr		
Platform Dimension (LxDxH)	380x300 mm		
Overall Dimension (LxDxH)	705x515x470 mm		
Weight	38 kg	38 kg	46 kg
Power	1600 W		
Power Supply	220V, 50Hz		
Catalog No.	9501283127	9501284127	9501285127

Shaking Water Bath LSWB-200 Series



Labocon Shaking Water Bath LSWB-200 Series is equipped with microprocessor controller for temperature and speed. It offers accurate heating up processes and simultaneous gentle mixing or vigorously shaking to agitate samples.

Features:

- * Digital display.
- * Stainless steel internal chamber.
- * Easy to clean.
- * Equipped with clamp.



Application:

Shaking Water Bath LSWB-200 Series from Labocon is used in biology, biochemistry, industry and pharmacology applications. Their scope of applications also comprises in incubation, fermentation, homogenization, chemical reactions, enzyme studies and tissue studies.

Specification:

Model	LSWB-201	LSWB-202	LSWB-203
Capacity	32 Litres		
Method	Orbital	Reciprocal	Orbital and Reciprocal
Temperature Range	RT+5°C to 65°C		
Speed Range	0-220 rpm	0-300 rpm	0-220 rpm; 0-300 rpm
Orbit Diameter	30 mm	20 mm	30 mm; 20 mm
Timer Range	0-999 hr		
Platform Dimension (LxDxH)	430x340 mm		
Overall Dimension (LxDxH)	705x455x440 mm		
Weight	38 kg	38 kg	46 kg
Power	400 W		
Power Supply	220V, 50Hz		
Catalog No.	9501286127	9501287127	9501288127



LABOCON SYSTEMS

3rd Floor, 207 Regent Street
London, W1B 3HH, UK
Email: info@labocon.com,
Web: www.labocon.com